


ASPIRATION BAG SYSTEMS



About Us

Founded in 2003 as Olgun Medical Construction and Trade LTD. ŞTİ, our company started to serve the medical sector with the sale of radiology consumables. Our company, which continues to sell other medical supplies, hospital equipment and various medical supplies in addition to the radiology consumables it specializes in, has adopted the basic principle of honesty, continuity, protecting customer interests and respecting human rights in service.

Our company, with its staff organized in the fields of purchasing, marketing, after-sales services and logistics support, prepares the ground for the improvement and perfection of the increasing service expectations of public and private health institutions, food and packaging sector in the field of gloves and consumables, and considers it as a basic duty to operate as an institution that has made its name in the health, food and packaging sector for years in order to provide quality products to its customers.

Our company has started its production activities by establishing a production facility in accordance with medical product standards for the production of Aspiration Bag Set and Canister product as the first stage, which has been continuing its R&D studies for about 3 years, and other consumables planned later.

Our vision is to increase the quality of life of people based on the power of teamwork and to be the leading company in the health sector, especially in our country and all over the world.

Our Mission; Olgun Medikal has a mission that aims to be a leader in the rapidly developing health, food and packaging sector with a long-term, permanent, customer satisfaction based, “professional service understanding and uncompromising quality” approach.

Our Values; Our VALUES, which represent the set of beliefs that make our company unique and the cultural characteristics that define the character of our company, are as follows; Customer Oriented, Reliable, Innovative, Responsible and Dedicated.

Our strength is our solidarity and team spirit.



As Olgun Medikal, in our activities in the medical sector;

- We aim to serve in the medical sector in accordance with national and international standards, regulations, and requirements,
- By keeping customer satisfaction at the highest level,
- By delivering products on time,
- By continuously improving the quality management system,
- By after-sales customer follow-up with our experienced and competent staff,
- By our sales and service activities in line with the latest technological developments,
- With the importance we attach to human health.

Standards

As Olgun Medical, a quality management system has been established in accordance with all standards required by the Medical Device Regulation (MDR EU 2017/745) directive published in May 2017, with certificates (TSE EN ISO 9001:2015 and EN ISO 13485:2016), in accordance with the international quality system standard, and our work continues accordingly.

Human Resources

Olgun Medikal, which aims to adhere to the principle of strong employees, strong company, has been continuing its way steadily since 2003 with its young and dynamic staff in a system where the performance levels of trained, self-educated individuals are monitored and the contributions that can be made to the processes are discussed.



DISPOSABLE ASPIRATION BAG SYSTEM

The Omsuction single vacuum bag system is designed for “high flow, high vacuum applications”. This system is the best in terms of practicality and safety in use. It is used both for operators and for hospitals and the environment. Ergonomic and robust construction ensures safe waste storage. The disposable suction bag system reduces the possibility of leaks and leakage directly to zero thanks to the strong welded structure of the bag and lid. It has successfully passed all relevant standard tests.





ASPIRATION BAG AND CANISTER JAR

Suction Aspiration Bag

- It aims to be indispensable for operations with its high suction power and superior material structure.
- With the plugging ports placed on the cover, it can be converted into safe waste without the need for extra accessories.
- It aims to minimize the risk of contamination in the environment based on the sealing power in the closed system.
- Designed in 1000 ml, 1200 ml, 2000 ml and 3000 ml capacities to meet all needs.

Suction Canister Jar

They are 1000 ml, 2000 ml, and 3000 ml capacity reservoirs that are connected to the vacuum regulator and placed inside the Disposable Aspiration Bag.

Thanks to its durable polycarbonate material structure, it is integrated with the bag and prevents waste from spreading to the outside environment.

It is cleaned in 121 degree sterilization devices.

The amount of liquid is measured with the volume indicators on it.

ASPIRATION BAG



OM0030



OM0031



OM0032

OM0030 Omsuction Aspiration Bag Set 1000 ml

OM0031 Omsuction Aspiration Bag Set 2000 ml

OM0032 Omsuction Suction Bag Set 3000 ml

ASPIRATION JAR



OM0014



OM0016




OM0015

Omsuction Aspiration Bag Jar 1000 ml **OM0014**

Omsuction Aspiration Bag Jar 2000 ml **OM0015**

Omsuction Aspiration Bag Jar 3000 ml **OM0016**

Aspiration Bag and Jar





ACCESSORIES

Accessoires



OM0006



OM0007

OM0007

Omsuction Tandem Elbow

OM0006

Omsuction Patient Elbow

OM0017

Omsuction Canister Elbow

OM0018


Omsuction T Connection Elbow



OM0017
OM0018



OM0020



OM0019



OM0025

OM0019 Omsuction Vacuum Directional Valve

OM0020 Omsuction Wall Mount

OM0025 Omsuction Rail Mounting






ASPIRATION BAG SET AND MOBILE TROLLEY

Suction Bag Set and Mobile Carrier Trolley

- Thanks to its ergonomic structure, it provides comfortable use in hospital conditions.
- Full capacity reaches approximately 12 L when connected with 4x3l jars.
- The jars can be fixed directly on the mobile trolley without the need for any other parts.
- When desired, the wheels of the carrier foot can be made brakeable.



**Maximum
Suction Power
performance**




**Secure Waste
Storage**



Suction

Aspiration
Bag Systems
Whole Set



CERTIFICATE



OLGUN MEDİKAL İNŞAAT VE TİCARET LİMİTED ŞİRKETİ

MAHMUTBEY MAHALLESİ ORDU CAD. NO: 29/11 BAĞCILAR - İSTANBUL - TÜRKİYE

Sales and Technical Service of Aspiration Bags and Cannisters; Sales of Operating Room Equipment and Medical Consumables; Sales, Installation and Technical Service of Radiology Imaging Medical Devices, Endoscopy Imaging Devices and Consumables

with a scope of

EN ISO 13485:2016

Has established a management system in accordance with international Medical Devices Quality Management System Standard

"Following elements of the standard are excluded"
"7.3" "7.5.5" "7.5.7" "7.5.9.2"

Certificate No	: M 11925
Initial Certification Date	: 26 May 2023
Certification Date	: 26 May 2023
Expiration Date	: 25 May 2026

General Manager



TÜRKAK BDS NO
YS-2084-4A7F

Last Modified : 26 May 2023 - R 00

Kiwa Belgelendirme Hizmetleri A.Ş.
İTOSB 9. Cadd No. 15 Tepeören Tuza
İstanbul / Turkey

Tel: + 90 216 593 25 75
Faks: + 90 216 593 25 74
info@kiwa.com.tr
www.kiwa.com.tr

Certificate is valid till expiration date,
subject to successful completion of
periodical surveillance audits.
Please contact above numbers for
detailed information.

CERTIFICATE



OLGUN MEDİKAL İNŞAAT VE TİCARET LİMİTED ŞİRKETİ

MAHMUTBEY MAHALLESİ ORDU CAD. NO: 29/11 BAĞCILAR - İSTANBUL - TÜRKİYE

Sales and Technical Service of Aspiration Bags and Cannisters; Sales of Operating Room Equipment and Medical Consumables; Sales, Installation and Technical Service of Radiology Imaging Medical Devices, Endoscopy Imaging Devices and Consumables

with a scope of

ISO 9001:2015

Has established a quality management system in accordance with international standard.

"Following elements of the standard are excluded"
"8.3"

Certificate No	: M 11924
Initial Certification Date	: 26 May 2023
Certification Date	: 26 May 2023
Expiration Date	: 25 May 2026

General Manager



TÜRKAK BDS NO
YS-ABE3-30A2

Kiwa Belgelendirme Hizmetleri A.Ş.
İTOSB 9. Cadd No. 15 Tepeören Tuza
İstanbul / Turkey

Tel: + 90 216 593 25 75
Faks: + 90 216 593 25 74
info@kiwa.com.tr
www.kiwa.com.tr

Certificate is valid till expiration date,
subject to successful completion of
periodical surveillance audits.
Please contact above numbers for
detailed information.

Last Modified - 26 May 2023 - R 00



OLGUN MEDİKAL İNŞ. ve TİC. LTD. ŞTİ.

Head Office:

Mahmutbey Mahallesi, Ordu Caddesi,
Perkotek Plaza No:29/1 Bağcılar/İstanbul-Türkiye

Factory:

Pelitli Mahallesi 4438 Sok 33/C Gebze / KOCAELİ

Tél: +90 212 438 46 71 info@olgunmedikal.com
Faks: +90 212 438 50 01 www.olgunmedikal.com