

CONSUMABLES



om **OLGUN**[®]
MEDİKAL



Founded in 2003 as Olgun Medikal our company began serving the medical sector for the purpose of radiology consumables sales.

Our company also continues to sell other medical equipment, hospital hardware and various medical supplies besides radiology consumables which it specializes in. It has adopted the basic principle of honesty in service, sustainability, protecting customer interests and respecting human rights.

Our company, with its organized personnel in the fields of purchasing, marketing, post-sales services and logistics support, provides a background in the fields of gloves and consumables for the improvement and perfection of the expanding service expectations of public and private health institutions, food and packaging industry and considers a fundamental duty to operate as an institution that made a name in the health, food and packaging sector for years in order to provide quality products to its customers.

Our vision;

To increase people's quality of life based on the power of teamwork, to become the leading company of the health sector worldwide, especially in our country.

Our mission;

As Olgun Medikal, we have a mission aiming to be a leader in its field with an understanding that "does not compromise on professional service context and quality" based on long-term, permanent, customer satisfaction in the health, food, and packaging industry.

Our values;

Our VALUES representing the sequence of beliefs that make our company unique and the cultural characteristics that define the character of our company are as follows; We are Customer Oriented, Reliable, Innovative, Responsible, and Dedicated.

Our strength is our solidarity and team spirit.

CONSUMABLES

MANUFACTURING

- ASPIRATION BAG SETS
- LABORATORY CONSUMABLES
- STERILIZATION REELS & POUCHES
- COLOSTOMY BAG & ADAPTER
- SUCTION CATHETERS
- YANKAUER SET

BRANDED

- MEDICAL GLOVES
- IV CANNULAS
- 3-WAY STOP COCK
- SPINAL NEEDLES
- ELECTROSURGICAL PENCIL
- ECG ELECTRODES
- IV INFUSION SET
- BLOOD TRANSFUSION SET
- URINE BAGS
- SUCTION CATHETERS
- FOLEY CATHETERS
- SURGICAL BLADES
- MEDICAL PAPERS



ASPIRATION BAG SETS

Omsuction Aspiration Bag Set

- Available size 1000mL, 2000mL, 3000mL with hose or w/o hose

Canister Jar

- Available sizes: 1000mL, 2000mL, 3000mL

Accessories

- Tandem elbow
- Patient elbow
- Canister elbow
- T Connection elbow
- Vacuum directional valve
- Wall bracket
- Rail bracket
- Mobile trolley








LABORATORY CONSUMABLES

WITH CONICAL BOTTOM TEST TUBE



TEST TUBE PS



TEST TUBE PP



EPPENDORF TUBE



INTEGRATED SYSTEM VACUUM
URINE SAMPLE CONTAINER



PNEUMATIC SYSTEM URINE
SAMPLE CONTAINER



BLUE PIPETTE TIP




PETRI DISH



YELLOW PIPETTE TIP



TUBE STOPPER FLANG TYPE



GAITA CONTAINER WITH CAP



URINE SAMPLE CONTAINER NON STERILE



GAITA CONTAINER WITH SPOON





STERILIZATION ROLL & POUCH



Sterilization Flat Reels

- Medical Paper provides maximum safety.
- All imprint are located outside the packaging area to prevent ink pigment migration to the product.
- Transparent , multilayer co-polymer film allows easy identification.
- Reinforced film to avoid tear during opening.
- Superior barrier and quality with 60gsm or 70gsm medical kraft paper.
- Indicators which are water based, non toxic and give accurate result for steam, EO Gas and Formaldehyde (FO).
- Strong seal-strength, visible sealing lines.
- Clean fibre-free opening.
- All materials comply with international standards.
- Opening direction marked.
- Conformity with EN 868-5.
- Available sizes: Any dimension between 5cm – 60cm can be made on request. On request, 100 m lenght can be produce as well.
- Sterilization Rolls can be produced with your trademark.



Sterilization Flat Pouches

- Medical Paper provides maximum safety.
- Pouches constructed from high grade medical kraft paper and multilayer transparent co-polymer film.
- Strong seal strength.
- Superior barrier and quality with 60gsm or 70gsm medical kraft paper.
- Indicators which are water based , non toxic and give accurate result for steam, EO Gas and Formaldehyde (FO)
- Strong seal-strength, visible sealing lines.
- Clean fibre-free opening.
- All materials comply with international standards.
- Opening direction marked.
- Conformity with EN 868-5.
- Apart from those standard sizes, as long as the width is one of the standards, there is no limit at lengths.
- Sterilization Flat Pouches are packed as 100 pieces wrapped in film according to requested dimension.
- Note: Sterilization Bags can be produced with your trademark.

COLOSTOMY BAG & ADAPTER



Designed for collecting waste from the colon (stoma) created on the abdominal wall due to bowel diseases (such as cancer, Crohn's disease, ulcerative colitis) or trauma.

- Consists of a skin-adhering base plate (adapter) and pouch system
- Prevents leakage and minimizes skin irritation
- Two-piece system providing secure and discreet ostomy protection
- Skin-friendly adhesive base plate with leak-proof sealing and odor control, designed for real-life use
- Base plate options for different stoma sizes (30, 45, 55, 70 mm)
- Base plate and pouch integrated with a locking system
- Pouch options available to meet different patient needs



SUCTION CATHETERS



- Sterile, single-use
- Rounded tip
- Female connector
- The color that makes the dimensions easy to identify
- Available Sizes: 6Ch / 8Ch / 10Ch / 12Ch / 14Ch / 16Ch / 18Ch

YANKAUER SET



- Kink-resistant PVC tubing
- Blue-colored connector
- Rigid and durable suction tip
- Cut-to-fit (cuttable) connection connector
- Large inner diameter to ensure unobstructed flow
- Four lateral holes





GLOVES



LATEX GLOVES POWDERED

Powdered Latex Examination Gloves

- Non-sterile
- Natural latex
- Single-use, powdered
- 100 gloves per box
- 20 boxes per carton
- 2,000 gloves per carton
- Available sizes: S, M, L



LATEX GLOVES POWDER FREE

Powder-Free Latex Examination Gloves

- Non-sterile
- Natural latex
- Single-use
- Powder-free
- 100 gloves per box
- 20 boxes per carton
- 2,000 gloves per carton
- Available sizes: S, M, L



NITRILE GLOVES POWDER FREE

OM Blue Nitrile Gloves are suitable for use in many industries, especially in healthcare institutions, workplaces, food, and beverage production with their thick feature. Not containing latex, they are preferred for latex-sensitive people.

- Nitrile
- Non-sterile
- Blue color
- Single-use
- Powder-free
- 100 gloves per box
- 20 boxes per carton
- 2,000 gloves per carton
- Available sizes: S, M, L



STERILE SURGICAL GLOVES

- Sterile (R), Single-use
- Medical paper on both sides of the unit packaging
- Powdered & Powder-Free
- 50 gloves per box
- 400 boxes per carton
- Available sizes: 6 / 6.5 / 7 / 7.5 / 8 / 8.5

IV CANNULAS



Mediflon™ / Super Life

IV Catheter with Port and Wing

It has all the specifications of a standard intravenous catheter.

Available sizes: 14G to 26G



Neocan™

IV Catheter for Babies and Newborns

Miniature IV catheter with small, thin, and smooth design especially for use in babies (children under 1 year) newborns (children under 1 month).

Size: 24G Port-free



Enfancie™

IV Catheter for Newborns

Miniature IV catheter with small, thin, and smooth design especially for use in newborns (children under 1 month).

Size: 26G Port-free



3-WAY STOPCOCK



- Sterile, single-use
- 360° rotating stopcock
- Flow direction indicated by arrows
- One male and two female Luer lock connectors
- Blue color
- Packing: 50 pieces per box / 1.000 pieces per carton

SPINAL NEEDLES



- Medispine Spinal Needle (Quincke Tip)
- Designed for spinal (lower Arachnoid) anesthesia.
- Available Sizes: 18G, 19G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G

ELECTROSURGICAL PENCIL



- Sterile, Single-use
- 3m cable
- 10 copper wires in each cable
- 3 pin plug
- Hand operated, yellow and blue
- 2 buttons

ECG ELECTRODES



ECG Electrode- Adult

- Single-use
- 50x55mm Medical grade foam
- With AgCl sensor
- Latex-free
- Excellent signal transmission
- high tack, short/long-term monitoring
- In special bags of 50 pieces

ECG Electrode- Pediatric

- Single-use
- 25,4x25,4mm Medical grade foam
- With AgCl sensor
- Latex-free
- Excellent signal transmission
- high tack, short/long-term monitoring
- In special bags of 50 pieces

STRESS TEST ECG ELECTRODE



Intended for use in applications where high tack adhesive is needed.

Recommended usage: ICU/CCU, Surgery, Stress, Holter, Emergency Cloth, Adult, Solid Gel, Round, 50mm
50pcs/pk, 40pk/cs

IV INFUSION SET



- Sterile, Single-use
- Small drip chamber
- Suitable for infusion bottles and bags
- Filter in the drip chamber
- Roller drop adjuster
- Covered airway
- Luer Lock
- 21G hypodermic needle

BLOOD TRANSFUSION SET



- Sterile, Single-use
- Non-toxic, pyrogen free
- Roller drop adjuster
- With blood filter
- With air flap
- 18G hypodermic needle




SUCTION CATHETERS



- Sterile, single-use
- Rounded tip
- Female connector
- The color that makes the dimensions easy to identify
- Available Sizes: 6Ch / 8Ch / 10Ch / 12Ch / 14Ch / 16Ch / 18Ch

FOLEY CATHETERS



- An Indwelling catheter used for short term urine drainage.
- Symmetrical large capacity balloon ensuring a straight tip for proper flow and good Sphincter action thus preventing leakage.
- Hard valve ensures easy inflation & deflation of balloon.
- Colour coded for instant size identification.
- Sterile, individually packed.
- Available in Two way and Three way Variants.
- 100% Silicone Coated
- Available sizes of Two Way: 6Fr / 8Fr / 10Fr / 12Fr / 14Fr / 16Fr / 18Fr / 20Fr / 22Fr / 24Fr
- Available sizes of Three Way: 14Fr / 16Fr / 18Fr / 20Fr / 22Fr / 24Fr

SURGICAL BLADES



- Sterile, Single-use
- Carbon steel
- Individually packed
- Boxes of 100 pcs
- Available Sizes: 10 – 11 – 12 – 15 – 18 – 20 – 21 – 22 – 24 – 25

URINE BAGS



Urine Bag - Adult

- Made of medical-grade soft PVC
- Non-toxic
- Equipped with an anti-reflux valve system
- Single-use
- Equipped with hanging holes
- Available Types: With T-tap, push pull valve, screw valve
- Volume: 2000mL



Urine Bag - Pediatric

- Single-Use
- Made of medical-grade soft PVC
- Available Volumes: 100mL and 200mL



Urometer

- Medical-grade PVC, ml graduated
- 2600mL urine bag, 500mL plastic chamber
- Hydrophobic bacterial filter, anti-reflux valve
- By-pass system, wide-bore and kink-resistant tubing
- Hanger with flexible fastening strap






MEDICAL PAPERS

UMUR Medical Papers

ECG, NST, STRESS TEST, DEFIBRILLATOR, EEG,
TYMPANOMETER, AUTOREFRACTOMETER, ULTRASOUND

umur



SONY Medical Papers

- UPP-110 S
- UPP-110 HG
- UPP-110 HD
- UPP-210 HD
- UPP-210 SE
- UPP-84 HG
- UPC-21 L (Colored)
- UPC-21 S (Colored)





OLGUN MEDİKAL İNŞAAT VE TİCARET LİMİTED ŞİRKETİ

Genel Merkez / Headquarter: Güzelyalı Mahallesi Sahil Bulvarı no:183/A Pendik/İstanbul
Fabrika Adres / Factory Address: Pelitli, Yeniçeşme Mevkii, Dostlar Cd. No:33, 41400 Gebze/Kocaeli

Tel : +90 216 519 29 77-79
Fax: +90 216 519 29 09

info@olgunmedikal.com
www.olgunmedikal.com